**bibag**On-line Dry Bicarbonate Concentrate





# bibag - the on-line dry bicarbonate concentrate for haemodialysis

for Fresenius Medical Care haemodialysis machines

#### The best solution

Use all the advantages, that bibag provides:

### Dry concentrate

- Small and lightweight
  - Easy handling
  - Small storage requirements
  - Easy transportation
- Chemically stable
  - Long shelf life

# On-Line dry concentrate

- Dissolved directly at the machine immediately before use
  - Prevention of microbiological growth
  - Safe dialysis treatment

# Unique additional advantages

- Flexible bag
  - Minimal waste volume
  - Comfortable to dispose of
  - Environmentally friendly
- Automatic Emptying Programme after treatment
  - No spilling or dripping during disconnection
  - No manual removal of remaining water
  - Most hygienic way to work



| Specifications       |              |                |
|----------------------|--------------|----------------|
| Packaging:           | units/carton | cartons/pallet |
| bi <i>b</i> ag 650 g | 16           | 56             |
| bi <i>b</i> ag 900 g | 12           | 56             |

Material bibag: PP/PE, Overwrap: PE, Shelf life: 3 years, CE-marked

| Order information |                         |  |
|-------------------|-------------------------|--|
| Weight/Bag        | Art. No.                |  |
| 650 g             | 508 992 1               |  |
| 900 g             | 508 991 1               |  |
| 650 g             | 506 078 1               |  |
| 900 g             | 506 080 1               |  |
|                   | 650 g<br>900 g<br>650 g |  |

\*bibag 5008 is the standard module for the 5008, 5008S and 4008S classic HD machines



